



Yeastern Biotech Co., Ltd

Product Use Limitation & Warranty

This product is intended to be used for life science research only. It has not been approved for drug or diagnostic purpose. YEASTERN's products should not be resold, modified for resale, or used to manufacture commercial products without written approval by YEASTERN. YEASTERN guarantees the performance of all products in the manner described in our protocol. The purchaser must determine the suitability of the product for its particular use. Should any product fail to perform satisfactorily due to any reason other than misuse, YEASTERN will replace it free of charge.

Ver. L0619

No part of these protocols may be reproduced in any form or by any mean, transmitted, or translated into a machine language without the permission of YEASTERN BIOTECH CO., LTD.

Address: 6F-3, 23 Lane 169, Kang Ning St.,
Shijr, Taipei, 22180 Taiwan.

Tel: +886-2-2695-3922 **Fax:** +886-2-2695-3979

Email: yeastern@yeastern.com.tw



O'in 1 DNA Polymerase Premix

Copyright© 2011 All rights reserved. Yeastern Biotech Co., Ltd.

Copyright© 2011 All rights reserved. Yeastern Biotech Co., Ltd.

Cat. No.
FYT201-100P
FYT202-100P

O'in 1 DNA Polymerase Premix

Volume: 5 ml/2.5 ml

Storage: -20°C

Description

O'in 1 DNA Polymerase Premix is an economical and ready-to-use premix, containing YEAtaq DNA Polymerase (# FYT001-500U), dNTP and all other reagents necessary for PCR, except DNA template and primers. It saves the time for preparing the master mix and reduces the risk of contamination from multiple pipetting steps.

Optimal PCR conditions, including template and primer concentrations and PCR program, should be determined experimentally by the investigator from case to case.

Component	FYT201-100P	FYT202-100P
O'in1 DNA Polymerase Premix	1X	2X

1x O'in1 DNA Polymerase Premix Contents:

10 mM KCl, 2 mM MgSO₄·7H₂O, 20 mM Tris-HCl (pH 8.8), 0.1 % Triton X-100, 10 mM (NH₄)₂SO₄, 0.1 mg/ml BSA, 0.2 mM dNTP mix, 50 U/ml YEAtaq DNA Polymerase, stabilizers

2x O'in1 DNA Polymerase Premix Contents:

Two folds concentration of all the above reagents

Procedure

A. Preparation of the PCR Master Mixes

1. Prepare a master mix according to Table 1 or 2.

Table 1. Reaction components when 1X O'in1 PCR Premix is used (FYT201-100P)

Component	Volume	Final conc.
1X O'in1 DNA Polymerase Premix	45 µl	1X
Forward Primer (10 µM)	1 µl	0.2 µM
Reverse Primer (10 µM)	1 µl	0.2 µM
Template DNA	0.5-3 µl	
ddH ₂ O	total to 50 µl	

Table 2. Reaction components when 2X O'in1 PCR Premix is used (FYT202-100P)

Component	Volume	Final conc.
2X O'in1 DNA Polymerase Premix	25 µl	1X
Forward Primer (10 µM)	1 µl	0.2 µM
Reverse Primer (10 µM)	1 µl	0.2 µM
Template DNA	0.5-10 µl	
ddH ₂ O	total to 50 µl	



2. Mix the master mixture thoroughly by pipetting up and down. Dispense the mixture into PCR tubes or plates.

B. Performing PCR

1. Program your instrument according to Figure 1.
2. Place the PCR tubes or PCR plates in the thermo cycler and start the cycling program.




Figure 1. O'in 1 PCR cycling conditions



Unit Definition

One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTPs into acid-insoluble material in 30 minutes at 72°C.

Quality Control

Nuclease activity is not detected after incubation 1µg lambda/HindIII DNA with 5 units O'in 1 DNA Polymerase in 50 µl reaction volume in supplied reaction buffer for 18 hours at 37°C.