

100 bp YEA Ladder DNA Marker IV

Ready-to-Use DNA Ladder

Recommended storage condition: -20 °C

Loading Volume: 2~3 μ L/well

Size Range: 100 bp~5 kb

Package size: 1 mL

The Yeastern Biotech's 100 bp YEA Ladder DNA Marker IV (FYD012-1ML) contains 15 discrete fragments ranging from 100 bp to 5000 bp. It is supplied in a ready-to-use format containing blue tracking dyes, and mixed at varied concentrations.

For best results, load 3 μ L of the marker per well. There are two fragments, 500 bp and 1200 bp, as the reference bands which are approximately 90 ng each for these two bands and 30 ng for other bands if load 3 μ L per well.

It allows sizing and concentration estimation of DNA fragments on agarose gels generated by PCR or restriction digestion.

E: yb@yb-biotech.com

T: +886-2-26953922


W: https://si.secda.info/yb_md/



Cat. NO. FYD012-1ML

100 bp YEA Ladder DNA Marker IV

Ready-to-Use DNA Ladder

	Size	ng/5 μ L
	5 kb	50
	3 kb	50
	2 kb	50
	1.5 kb	50
	1.2 kb	150
	1 kb	50
	900 bp	50
	800 bp	50
	700 bp	50
	600 bp	50
	500 bp	150
	400 bp	50
	300 bp	50
	200 bp	50
	100 bp	50

YB Yeastern
Biotech
Biotech for Better Life

Picture:

1.5% Agarose,
5 μ L/lane
1x TAE buffer